

Cryz J- Wb, Elisa - Q08257

مصنوع: CM0019456



مقدم

Cryz J- is a monoclonal antibody (mAb) derived from a transgenic mouse (KO) that produces human IgG1 antibodies. It is used for the detection of Cryz J- mRNA, protein, and antigen in various samples. The antibody is produced by a transgenic mouse (KO) that produces human IgG1 antibodies. It is used for the detection of Cryz J- mRNA, protein, and antigen in various samples. The antibody is produced by a transgenic mouse (KO) that produces human IgG1 antibodies. It is used for the detection of Cryz J- mRNA, protein, and antigen in various samples.

ذات صلة

اسم المنتج	معلومات
مستضد	Cryz J- Rabbit pAb
مستضد	Cryz J-; YZ
مستضد	مستضد
مستضد	1429
مستضد	1429
مستضد	NP_001123514.1 (Cryz J- Rabbit pAb) 1-329 amino acid sequence.
مستضد	مستضد
مستضد	مستضد
مستضد	WB, ELISA
مستضد	مستضد
مستضد	Q08257
مستضد	35kDa
مستضد	مستضد