

ApC Anti-Human Cd335/Nkp46 Monoclonal Antibody For Flow Cytometry - 076036

رصدن عل مقرر: CM0005913



مقدم

APC-conjugated rabbit monoclonal antibody targeting human CD335/NKp46 (SWISS: 076036) specifically for flow cytometry. Reacts with human NK cells, useful for tumor cell lysis studies.

ديزل فرعا

راي عمل	همل
جتنم الم س	APC Rabbit anti-Human CD335/NKp46 mAb
راعستم مس/اتاطحال م	LY94; CD335; NKp46; NK-p46
عاون ال	Human
نجل الفرع (ناسن)	9437
نجل الفرع	9437
دستم ال	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-254 of human CD335/NKp46 (NP_004820.2)
ردصم ال	Rabbit
هفل	Monoclonal Antibodies
قبطط ال	FC
عطاقم ال لعافت ال	Human
SWISS	076036
نيتورب ال نرو	21kDa/22kDa/23kDa/32kDa/34kDa
نشل ال	Ice bag