

Mouse Anti-Mouse/Rat Cd90.1 Monoclonal Antibody For Fc - P01831

رخصن عل مقرر: CM0008298

مقدم

-conjugated mouse anti-mouse/rat CD90.1 (Thy-1) monoclonal antibody for flow cytometry. Detects the GPI-anchored membrane glycoprotein involved in synaptogenesis and cell-cell interactions. Suitable for neuronal and immunological research in mouse and rat.

دی زمل فرعا

Parameter	Value
Product Name	Mouse anti-Mouse/Rat CD90.1 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
GeneID (Human)	7070
GeneID	21838
Immunogen	Recombinant protein of mouse/rat CD90.1
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P01831
Shipping	Ice bag