

Mouse Anti-Rabbit IgM (Fc) Monoclonal Antibody For Flow Cytometry - P03988

رصدن عل مقرر: CM0006983

مقدم

Mouse monoclonal antibody (conjugate) specific for rabbit IgM Fc region. Recombinant immunogen (aa 107-458, P03988). Validated for flow cytometry. Ideal for diagnostic assay development and research.

دیزمل فرع

Parameter	Value
Product Name	Mouse Anti-Rabbit IgM (Fc) mAb
Species	Rabbit
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 107-458 of rabbit IgM (Fc) (P03988)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P03988
Protein Weight	52kDa
Shipping	Ice bag