


Pe Mouse Anti-Rabbit IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01870

رصدن عل مقرر: CM0006472



مقدم

PE-conjugated mouse monoclonal antibody targeting the Fc region of rabbit IgG. Designed for flow cytometry detection of rabbit immunoglobulins in immunological research and diagnostic assays.

دیزمل فرعا

Parameter	Value
Product Name	PE Mouse anti-Rabbit IgG (Fc) mAb
Species	Rabbit
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 96-323 of rabbit IgG (Fc) (P01870)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P01870
Protein Weight	35kDa
Shipping	Ice bag