

Rabbit Anti-Human Egfr Mab For Fc - P00533

رصدن عل مقو: CM0008412



ةم دقم

Rabbit anti-Human EGFR monoclonal antibody validated for flow cytometry in human samples. Targets extracellular domain (aa 335-525), enabling receptor signaling studies and oncology research applications.

دي زم ل فرع

لم اعلم ل	هم ي قول
ج تنم ل م سا	Rabbit anti-Human EGFR mAb
م س/ا ط ا ح ال م راع س ت م	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2
ع و ن ل	ناس ن
GeneID (ناس ن)	1956
د ص ت م ل	يرش ب ل EGFR م 335-525 ه ن ي م أ ل ص ا م أ ل ع م ق ف ا و ت ي ل س ل س ت ي ل ع ي و ت ح ي ب ي ك ر ت ل ا د ا ع م ي ج ا م د ن ا ن ي ت و ر ب ا ب ي ك ر ت ل ا د ا ع م ن ي ت و ر ب (NP_005219.2).
ر د ص م ل	ب ن ر أ
ة ئ ف ل	ة ل ي س ن ل ا ه ي د ا ح ا ه ا ص م م ا س ج أ
ق ي ب ط ت ل	FC
ع ط ا ق ت م ل ا ع ا ف ت ل	ناس ن
SWISS	P00533
ن ي ت و ر ب ل ن ز و	44kDa/69kDa/77kDa/134kDa
ن ح ش ل	ج ل ت س ي ك