


Rabbit Anti-Human Vcam-1/Cd106 Mab For Fc - P19320

رصدن عل مقرر: CM0004837



مقدم

Rabbit anti-Human VCAM-1/CD106 monoclonal antibody for flow cytometry (FC). Validated for human samples; targets vascular cell adhesion protein 1 involved in immune surveillance and inflammation.

دیزمل فرع

لماعمل	مقوال
جتنملا مسا	Rabbit anti-Human VCAM-1/CD106 mAb
راعسم مس/اتاطحال م	CD106; INCAM-100
عونلا	Human
نجلال فرعم	7412
دصتسملا	Recombinant fusion protein containing a sequence corresponding to amino acids 25-698 of human VCAM-1/CD106 (NP_001069.1).
ردصملا	Rabbit
مفل	Monoclonal Antibodies
قبطتلا	FC
عطاقملا لعاقتلا	Human
SWISS	P19320
نرتوربلا نزو	71kDa/74kDa/81kDa
نحشلا	Ice bag