


Rabbit Anti-Mouse Cd274/Pd-L1 Monoclonal Antibody () For Flow Cytometry - Q9Ep73

رصنعلا مقر : CM0008928



مقدم

-conjugated rabbit monoclonal antibody against mouse CD274/PD-L1, validated for flow cytometry. Ideal for immune checkpoint studies and tumor immunology research.

ديزملل فرعلا

Parameter	Value
Product Name	Rabbit anti-Mouse CD274/PD-L1 mAb
Remarks/Alias	B7h1; Pdl1; Pdcd1l1; Pdcd1lg1; A530045L16Rik
Species	Mouse
GeneID (human)	29126
GeneID	60533
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9EP73
Protein Weight	32kDa
Shipping	Ice bag