


Yes1 Rabbit Pab For Wb, If/Icc, Elisa - P07947

رصدنعل مقرر: CM0014380



مقدم

YES1 Rabbit Polyclonal Antibody (pAb) validated for WB, IF/ICC, ELISA. Recognizes Human and Mouse Tyrosine-protein kinase Yes (c-yes), a regulator of cell growth, adhesion, and migration. Ideal for cancer research and signal transduction studies.

دیزمل فرعا

لماعمل	هميقل
جتنملا مس	YES1 Rabbit pAb
راغتسم مس/اتاطحال م	Yes; c-yes; HsT441; P61-YES; YES1
عاونالا	Human
GeneID (human)	7525
GeneID	7525
دصتسملا	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human YES1 (NP_005424.1).
ردصملا	Rabbit
ةفلا	Polyclonal Antibodies
قبططالا	WB, IF/ICC, ELISA
عطاقتملا لغافتلا	Human, Mouse
SWISS	P07947
نيتوربالا نرو	61kDa
نحشلا	Ice bag